

Inspired Care

For your Bi-level/CPAP Patients

RT040 Full Face Mask

Respiratory

HUMIDIFICATION

Fisher & Paykel
HEALTHCARE

RT040 Hospital Single Use Full Face Mask

The easy-fit mask solution with your patient's comfort in mind

Unique Headgear System

- Easy to secure
- Cushioned forehead rest

RT040S – Small
RT040M – Medium
RT040L – Large

FlexiFit™ Technology

- Designed to overcome nasal bridge breakdown

Exhalation Vents

- Minimal deadspace
- Quiet operation

Built-in Ports

- Pressure
- Oxygen

Dual Swivel Elbow

- Aids patient movement

Anti-Asphyxiation / Non-rebreathing Valve


- Prevents re-breathing
- Opens if ventilator fails

Free Motion Glider™ Strap

- Quick release
- Single motion refit
- Provides freedom of movement

Under Chin Seal

- Comfortable
- Soft & flexible
- Stable
- Silicone



Achieve optimal performance with
Fisher & Paykel Healthcare's MR850 Humidification System

INTERNATIONAL P O Box 14 348, Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com
AUSTRALIA Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **AUSTRIA** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **BENELUX** Tel: +31 78 644 0924 Fax: +31 78 644 0914
CHINA Tel: +86 20 3205 3486 Fax: +86 20 3205 2132 **FRANCE** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **GERMANY** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99
INDIA Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **IRISH REPUBLIC** Tel: 1800 409 011 **ITALY** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709
JAPAN Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **MIDDLE EAST** Tel: +9 71 4 360 0849 Fax: +9 71 4 368 8110 **SPAIN** Tel: +34 902 013 346 Fax: +34 902 013 379
SWITZERLAND Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Tel: +44 1628 626 136 Fax: +44 1628 626 146 **Nthn Ireland** Tel: 0800 132 189
USA Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

Respiratory

HUMIDIFICATION

Fisher & Paykel
HEALTHCARE